LESTBOLL Bath Analyzer 20

Area of application for the ZESTRON® Bath Analyzer 20




In order to achieve consistent cleaning results, the user has to ensure that the actual concentration of the cleaning agent stays within the recommended application range. Various factors such as "drag-out", dilution and evaporation can potentially influence the actual concentration. Regular measurements of all process related parameters are therefore highly recommended. The ZESTRON® Bath Analyzer 20 is a simple, easy to use method which provides reliable and accurate results for fresh as well as for contaminated cleaning baths.

Simple measurement in only a few steps:


① Sampling / Decanting







Shaking / Evaluating



Advantages of the ZESTRON® Bath Analyzer 20:

- Very precise concentration measurement:
 - even for increasing contamination levels
 - for fresh as well as for contaminated cleaning baths
- Simple and fast procedure
- Mobile, for use anywhere
- Requires little space with minimal investment

Specifically developed for the following cleaning agents:

- VIGON® SC 200
- VIGON® N 600
- VIGON® N 640
- VIGON® SC 202 • VIGON® SC 210
- VIGON® PE 180
- VIGON® TC 150
- VIGON® PE 200
- VIGON® PE 215N
- VIGON® PM 105
- VIGON® 1000 CR
- HYDRON® SE 220
- HYDRON® SE 230A
- HYDRON® SC 300
- ZESTRON® VD 200

The ZESTRON® Bath Analyzer 20 Kit includes:



- 1) PE bottle with injection nozzle
- 2) 2 pairs of vinyl gloves
- 3) Digital thermometer (°C/°F)
- 4) Sampling beaker
- 5) Glass cylinder with marking
- 6) 6 x 100 ml ZESTRON® Bath Analyzer 20 test solution
- 7) User manual, including correlation charts (not illustrated)

PR00001735.17.DOC Release Date: 30.04.2018

South Asia